

Veterinary clinic equipment

Hospital cages, infrared therapy cages, operating and diagnosis tables, bathing tanks, grooming tables and dryer rooms. Every model with its finished dimensions.

Items listed	Groups	Packaging
22	4	Own brand or private label

Cages & housing · 4

CAT. NO.	PRODUCT	MATERIAL	DIMENSIONS / CONFIGURATION
AHC-PHC03	Pet hospital cage, four compartment	Stainless steel	1220 × 752 × 1560 mm
AHC-PHC04	Pet hospital cage, four compartment, tall	Stainless steel	1220 × 700 × 1570 mm
AHC-PHC05	Pet hospital cage, six compartment	304 stainless steel	1500 × 700 × 1570 mm
AHC-PHC06	Pet hospital cage, six compartment, deep	304 stainless steel	1500 × 605 × 1700 mm

Grooming · 4

CAT. NO.	PRODUCT	MATERIAL	DIMENSIONS / CONFIGURATION
AHC-PGT03	Electric grooming table, oval	Stainless steel	1180 × 600 × 600–900 mm
AHC-PGT04	Electric grooming table, rectangular	Stainless steel	1200 × 600 × 685–990 mm
AHC-PGT02	Electric grooming table, round	Stainless steel	600 mm diameter, 680–1080 mm
AHC-PGT01	Electric grooming table, scissor lift	Stainless steel	1250 × 650 × 250–950 mm

Recovery & therapy · 4

CAT. NO.	PRODUCT	MATERIAL	DIMENSIONS / CONFIGURATION
AHC-PDR03	Automatic pet dryer room, large, twin door	Powder-coated steel	1200 × 760 × 1190 mm
AHC-PDR02	Automatic pet dryer room, medium	Powder-coated steel	810 × 730 × 990 mm
AHC-PTC02	Infrared therapy cage, four compartment	Stainless steel	1220 × 700 × 1830 mm
AHC-PTC01	Infrared therapy cage, oxygen capable	Stainless steel	1220 × 700 × 1830 mm

Tables & wet rooms · 10

CAT. NO.	PRODUCT	MATERIAL	DIMENSIONS / CONFIGURATION
AHC-POT02	Electric operating table, column, tilting	Stainless steel	1610 × 720 × 785–1090 mm, tilt 0–15°
AHC-POT01	Electric operating table, scissor lift	Stainless steel	1300 × 600 × 450–1050 mm
AHC-PDT04	Pet diagnosis table with scale	Stainless steel	1220 × 600 × 750–1050 mm
AHC-PDT02	Pet diagnosis table, column	Stainless steel	1195 × 600 × 750–1055 mm

Medicare Vet · Veterinary clinic equipment

Specifications subject to change · quotation on request

CAT. NO.	PRODUCT	MATERIAL	DIMENSIONS / CONFIGURATION
AHC-PDT05	Pet diagnosis table, pedestal	Stainless steel	1100 × 600 × 800 mm
AHC-PDT03	Pet diagnosis table, two tier	Stainless steel	1100 × 550 × 645 mm
AHC-PDT01	Pet disposal table with sink and drawers	Stainless steel	1130 × 550 × 840 mm
AHC-PWT03	Stainless steel bathing tank with ramp	Stainless steel	1300 × 743 × 1290 mm
AHC-PWT05	Stainless steel bathing tank, electric lift	Stainless steel	1298 × 740 × 1150–1830 mm
AHC-PLT01	Stainless steel luxury trolley	Stainless steel	896 × 493 × 856 mm

Availability and registration requirements differ by country. In the United States, veterinary tests for infectious disease are regulated as veterinary biologics by USDA APHIS and require a licensed importer; non-infectious veterinary tests fall under FDA CVM. Certificates are issued by the manufacturing site and supplied with the quotation. Registration and import licensing remain with the importer of record.