

Isothermal nucleic acid test cards

Chairside molecular testing on a palm-sized dock. Cold-chain-free reagent, up to six pathogens per card, result in under forty minutes.

Items listed	Groups	Packaging
17	4	Own brand or private label

Gastrointestinal · 5

CAT. NO.	PRODUCT	SPECIES	SPECIMEN	PACK
ISO-FPV-CPV	FPV / CPV Nucleic Acid Test Card FPV/CPV	Dog & cat	Faeces / rectal swab	Single-use card
ISO-FCoV	Feline Coronavirus Nucleic Acid Test Card FCoV	Cat	Faeces / rectal swab	Single-use card
ISO-GIA-TRF-PV-CoV	Giardia / T. foetus / Parvovirus / Coronavirus Nucleic Acid Test Card GIA/TRF/PV/CoV	Dog & cat	Faeces / rectal swab	Single-use card
ISO-GIA-TRF	Giardia / Tritrichomonas foetus Nucleic Acid Test Card GIA/TRF	Dog & cat	Faeces / rectal swab	Single-use card
ISO-PV-CoV	Parvovirus / Coronavirus Nucleic Acid Test Card PV/CoV	Dog & cat	Faeces / rectal swab	Single-use card

Notifiable disease · 1

CAT. NO.	PRODUCT	SPECIES	SPECIMEN	PACK
ISO-ASFV	African Swine Fever Virus Nucleic Acid Test Card ASFV	Swine	Whole blood / serum	Single-use card

Respiratory · 4

CAT. NO.	PRODUCT	SPECIES	SPECIMEN	PACK
ISO-CDV-CPIV-CAV-2-Bb-CIV-Myc	CDV / CPIV / CAV-2 / Bb / CIV / Mycoplasma Nucleic Acid Test Card CDV/CPIV/CAV-2/Bb/CIV/Myc	Dog	Eye, nose or throat swab	Single-use card
ISO-CDV	Canine Distemper Virus Nucleic Acid Test Card CDV	Dog	Eye, nose or throat swab	Single-use card
ISO-FCV-FHV-1-Mf-Cf-Bb	FCV / FHV-1 / M. felis / C. felis / Bb Nucleic Acid Test Card FCV/FHV-1/Mf/Cf/Bb	Cat	Eye, nose or throat swab	Single-use card
ISO-FCV	Feline Calicivirus Nucleic Acid Test Card FCV	Cat	Eye, nose or throat swab	Single-use card

Vector-borne & parasitic · 7

CAT. NO.	PRODUCT	SPECIES	SPECIMEN	PACK
ISO-ANA-EHR	Anaplasma / Ehrlichia Nucleic Acid Test Card ANA/EHR	Dog	Whole blood / serum	Single-use card
ISO-BG-BC	Babesia gibsoni / canis Nucleic Acid Test Card BG/BC	Dog	Whole blood / serum	Single-use card
ISO-DI-DR	Dirofilaria immitis / Dirofilaria repens Nucleic Acid Test Card DI/DR	Dog	Whole blood / serum	Single-use card
ISO-FeLV-FIV	FeLV / FIV Nucleic Acid Test Card FeLV/FIV	Cat	Whole blood / serum	Single-use card
ISO-fHM	Feline Haemotropic Mycoplasma Nucleic Acid Test Card fHM	Cat	Whole blood / serum	Single-use card
ISO-LSH	Leishmania spp. Nucleic Acid Test Card LSH	Dog	Whole blood / serum	Single-use card
ISO-TOXO	Toxoplasma gondii Nucleic Acid Test Card TOXO	Dog & cat	Whole blood / serum	Single-use card

Availability and registration requirements differ by country. In the United States, veterinary tests for infectious disease are regulated as veterinary biologics by USDA APHIS and require a licensed importer; non-infectious veterinary tests fall under FDA CVM. Certificates are issued by the manufacturing site and supplied with the quotation. Registration and import licensing remain with the importer of record.